

StoichiometryFitter: An Interactive WebApp for Mineral Analysis

Zack Gainsforth¹, Michael Ofengenden^{1*}, and Roger Yu^{1*}

¹ Space Sciences Laboratory, University of California Berkeley, 7 Gauss Way CA, USA * These authors contributed equally.

DOI: [10.xxxxxx/draft](https://doi.org/10.xxxxxx/draft)

Software

- [Review](#)
- [Repository](#)
- [Archive](#)

Editor: [Open Journals](#)

Reviewers:

- [@openjournals](#)

Submitted: 01 January 1970

Published: unpublished

License

Authors of papers retain copyright and release the work under a Creative Commons Attribution 4.0 International License ([CC BY 4.0](#)).

Summary

This project is based off a GUI we had created in 2016, for fitting solids and analyzing mineral phases. We had developed a software tool for the analysis of multi-phase quantification of materials. One would provide the quantification of elements as input in the form of atomic %, elemental weight %, oxide weight % (with cation oxidation states), or as raw counts of TEM measurements. When using raw counts, one could also use the GUI to apply a set of k-factors, a thickness correction, and a detector correction. More of the use cases of the 2016 version are listed in Gainsforth's Previous paper¹. In 2022, we worked on updating Stoichiometry Fitter to the current version, including updating to the latest python version, providing more capabilities, adding more phases, and integrating QMin² a mineral recognition webapp. Then by converting Stoichiometry Fitter from a GUI to web app, our hope is that this tool will be more readily available to the public and will help make the research of geologists, mineralogists, and chemists alike, more efficient and easy when exploring new minerals.

Statement of need

Stoichiometry Fitter allows a researcher to convert elemental quantification to meaningful conclusions about mineral phases through python scripts. Researchers will then be able to reproduce results on multiple inputs using the same python scripts. The program produces standard outputs for various minerals including Ternary diagrams, diagrams for MORE EX-AMPLES, allowing and improving reproducibility for Geologists and researchers on different elemental compositions. Stoichiometry Fitter also allows the user to remove common artifacts from these quantifications simplifying the understanding process and easing the procedure of mineral analysis. Stoichiometry Fitter is an easy to use program that provides researchers with detailed analyses of a specified mineral. In the case where the researcher is unsure of the mineral that is composed by inputs of element quantifications, our team has integrated QMin into Stoichiometry Fitter. This allows for researchers to identify and classify new or unknown minerals.

State of Field

Stoichiometry Fitter allows for the process of analyzing minerals from element quantifications in a much more efficient scale than other programs researchers have used previously for the same kind of analysis. In field of research, many use spreadsheets or excel to write csv files or analyze elemental composition from their inputs. This can be slow and quite inefficient. Stoichiometry Fitter resolves this issue by presenting a scalable program that can analyze element quantifications from a direct input or a large text file. The program can then output multi-dimensional analysis. Along with the visuals and other integrated features of

40 Stoichiometry Fitter, it allows for a faster and more thorough than many other computational
41 spreadsheet programs.

42 Math & Figures

StoichiometryFitter

Element	Counts	Charge
H	0	1
He	0	0
Li	0	1
Be	0	2
B	0	3
C	0	4
N	0	-3
O	4	-2
F	0	-1
Ne	0	0
Na	0	1
Mg	2	2
Al	0	3
Si	1	4
P	0	5
S	0	-2
Cl	0	-1

Phase Analysis:
☒ Olivine

Arbitrary Absorption:
☒ Titan Detector Window

TEM Thickness Correction
☒ 0 g/cm³ * nm 18 takeoff in deg

Apply k-factor for
☒ Titan 60 keV

Oxygen by Stoichiometry?
☒ Silicates

Input Type:
☒ Counts
☐ At %
☐ Wt %
☐ Ox Wt % %

Go!

Input data: Element Counts O 4.000 Mg 2.000 Si 1.000 Total: 7.000
 Arbitrary Absorption Correction used: Titan Detector Window
 k-factors used: Titan 60 keV Oxygen determined by stoichiometry Quantification results: Element At% Wt% Ox Wt %
 Valence k-factor O 56.772 45.186 0.000 0.000 1.179 Mg 29.683 35.890 59.515 2.000 0.932 Si 13.544 18.924 40.485 4.000 1.000
 --- Olivine Analysis --- Mg/(Mg+Fe) = 1.000 Mg/(M1+M2) = 1.000
 Cations per 4 oxygens: Element: # Mg: 2.091 Si: 0.954
 Total Cats: 3.046

Result:

43 and referenced from text using ??.

45 Chemistry Calculations:

46 Citations

47 pandoc -citeproc paper.bib

48 References

- 49 ■ name: Gainsforth, Z. "Stoichiometry Fitter, a GUI for Fitting Solid Solutions
50 and Analyzing Mineral Phases: Microscopy and Microanalysis." Cambridge Core,
51 Cambridge University Press, 25 July 2016, <https://www.cambridge.org/core/journals/microscopy-and-microanalysis/article/stoichiometry-fitter-a-gui-for-fitting-solid-solutions-and-analyzing-mineral-phases/84C98F55BAE56F9BAED27F1681B661C9>.
52 index: 1
- 53 ■ name: Silva, Guilherme Ferreira, et al. "Qmin: A Machine Learning-Based
54 Application for Mineral Chemistry Data Processing and Analysis." 2021,
55 <https://doi.org/10.21203/rs.3.rs-629516/v1>. index: 2